

PROGRAMME 50B: DISINFECTION BY-PRODUCTS IN HIGHLY MINERALISED MINERAL WATERS

150 € excl. VAT – total amount for 1 test (excluding transport costs)



New in 2023



Need to test another method, evaluate your staff?

Order additional test samples (parcel in its entirety): 80 € excl. VAT (excluding transport costs)

Parameters to analyse	Bot	tle	Number of	
	Volume	Number	measurements per parameter and per bottle	
23M50B.1 – Highly mineralised mineral water - sent in February 2023 - Refrigerated parcel				
BrO ₃ -, ClO ₂ -, ClO ₃ -	250 mL	2	2	

PARTICULARITIES

Test conducted for a minimum of 10 participants.

Recommended period to start the sample treatment (PRDT):			
time interval during which the quality of test materials is optimal (in number of days)			
BrO ₃ -, ClO ₂ -, ClO ₃ -	D ₀ +3		

D₀: Day the samples are sent to all the participants (for most proficiency tests, Tuesday)